# BiliCare<sup>™</sup>

# Neonatal Noninvasive Bilirubin Meter

#### **The BiliCare<sup>™</sup> Transcutaneous Bilirubin Meter**

The BiliCare System measures the level of bilirubin using a patented technology to provide accurate results. It's technology is based on spectroscopy.

The System measures the level of transcutaneous bilirubin in the newborn by transmitting light at di° erent wavelengths through the outer ear. The amount of light absorbed by the bilirubin is calculated according to a customized algorithm.

- Less sensitive to motion artifacts
- Not sensitive to di<sup>o</sup> erences in pressure applied or angle of device

### The LED advantage

LEDs do not deteriorate over time. No shift in power or in wavelength means no need for routine calibration. Optional calibration check allows user to conÿrm calibration when desired. The calibration check process is quick and easy & provides immediate results.

### **Communication features**

BiliCare includes communication features that may connect it to a PC and/or to a HIS via HL-7

### **BiliCare Viewer**

Unique S/W that displays the results in tables and graphs. These may be printed and exported.

Search by Patient ID1:		0		20	0			
Patient Mile Selfuer	(market)	States Dates	同間	I ALLAND	C. C. Lawrence	1000 B		1
8/196 321	Patient Big:	28-Dec-13 18:00	1.0 #9/0.	a state Percentite	BRCare 2/91	348279		
Area ID.	Parisetti KU2:	30-Dec-14 17:00	6.7 mg/di	23-th-3400 Percentile	Addaw 5.76	548249		1.00
fulined IDIa Canada		Barth Dates . (* Dec 13 at 22					-	
Notice 30.	Patriet 810	23-Itec-19 18-80	11.3 mg/dl.	<b>Falls-Shift Personality</b>	BRCas LN	548219	5	-
Artes dit:	Pytient Kill:	30 Dec 15 22 18	7.2 mg/dL	Sigh JSth Percantile	BRC81 1-70	241275	E - 2	
Area ID:	Patient Eld:	\$9-Dec-25 18:85	2.5 eg/d.	25th-50th Percentile	\$8Cev 3.76	548219	Non	
Paddent Bills Diere	1000000000	Barth Dates 47 Dec 25 10.00	m 20		10000	233.58	1	
Birte ID:	Putners \$10:	31-0ec-19.00:00	6.3 mg/d.	B-35th Proceeding	BRCars 5/8:	345249	BillCare	
Auros stor	Pytent Kig:	20 Dec 15 03:00	11.0 mg/d.	75th \$5th Percentile	86Care 1/8	243218	- B	
Aires 80;	Parlwort KU2:	39-Dec-15 13:00	13.0 mg/dl.	In the local division of the local divisiono	BRCare S.N.	348299		1
Byles ID:	Pytierii 810.	88-Dec; 12 18:00	8.2 mg/d.	75th 85th Parcantlin	BiCev 1.76	141229	-	
Printed Bills Bartitle	1	HIGH DATE: LITTLE-TS LITTLE	12	THE R. LEWIS CO.				12
8/9e XV:	Patient (U2)	15-Dec 15-15-01	3.5 mg/dL	a state Percentile	10Care 2/10	246279		-
Nortes 80;	Patient KID:	34-Dec115-15-18	6.9 mg/d.	State 25th Percentile	MCan 5.76	548276		
Arise ID-	Pulieri IDD	\$3-Dec-13 33-24	0.3 mg/d.	S-1555 Percentile	BáCare S.N.	141276		



# www.mennenmedical.com

Israel: Mennen Medical Ltd. 6 Ha-Kishon Street, Yavne. P.O. Box 102, Rehovot 7610002, Israel. Phone: +972-8-9323333. **Fax:** +972-8-9328510



Defining tomorrow, today.

United States: Mennen Medical Corp. 290 Andrews Road, Feasterville-Trevose, PA 19053-3480, USA. **Phone:** +215 259-1020. Fax: +215 357-2010

Contact Information Ruchi Tiwari Genworks Health Pvt Ltd Unit 522-524, 5th Floor, Gamma Block Sigma Tech Park, Ramagondanahalli Varthur Road, Bangalore – 560066 Ph no. +91 96204 06263 ruchi.tiwari@genworkshealth.com

United Kingdom: Charter Kontron Ltd. Unit 18 Avant Business Center 21 Denbigh Road Milton Keynes MK1 1DT, UK. Phone: +44 1908 646070. Fax: +44 1908 646030



# BiliCare<sup>™</sup>

# Neonatal Noninvasive Bilirubin Meter







Defining tomorrow, today

# Accurately and conveniently measures and monitors the level of bilirubin to assess the risk of neonatal hyperbilirubinemia.



## THE NEXT GENERATION TRANSCUTANEOUS BILIRUBIN METER

#### Simple for the User & Gentle on the Baby

- Reduces the need for Total Blood Bilirubin and the number of sticks in the baby's heel, by providing an accurate bedside measurement
- Only requires a single measurement per test

#### LED Transmission Technology Reduces User Variances & Maintenance Requirements

 Less sensitive to motion artifacts & measurement differences based on user technique • LEDs do not deteriorate over time eliminating the requirement for routine device calibration

## ADVANCED TECHNOLOGY FOR ACCURATE MEASUREMENTS

#### Speed and convenience

- Reduces number of blood tests via heel puncture
- Can use just one measurement or average 2-3 measurements
- · Can be used with or without entering baby's ID
- Large & clear display for easy use
- Resistive touch screen can be used with or without gloves
- Menu-driven user interface minimizes training needs, facilitating use by multiple caregivers
- Memory up to 40 measurements saves time in transcribing and comparing results
- Option to enter both patient and user identification facilitates hospital audits and quality assurance
- Barcode scanner for quick and accurate entry of caregiver and baby identification

#### Accuracy

- The BiliCare System was validated in clinical studies to produce results comparable to blood tests
- Repeatability tests confirm reliable results

#### Low cost of ownership

- The BiliCare System is capable of working with or without disposable tip cover
- · LEDs do not require routine calibration, minimizing maintenance and cost associated with service
- · Long-lasting, rechargeable battery with cradle allows extended use and ease of recharging



### **MEASURING WITH BILICARE**

Neonatal jaundice is a common condition, but when not carefully monitored, hyperbilirubinemia can lead to kernicterus, a devastating condition. Early measurement of bilirubin levels, leading to early intervention is recommended but the time to take a blood sample and the trauma to the baby made universal screening di cult . Now, the BiliCare System provides a cost and time e° ective way to monitor bilirubin levels without unnecessary heel sticks.

#### Ease of Use - just follow the simple prompts for quick measurement where and when you need it:

Simply press measure on the home screen to initiate the measurement process

Position the BiliCare clip on the baby's ear to take your measurement

You are now ready to review your results or take additional measurements per your hospital protocol and device set-up





