




MAGLUMI® X3

Fully-auto Chemiluminescence Immunoassay (CLIA) System



Save Your Space without Compromise

Powerful Performance and Compact Design Compatible with Small and Medium-sized Hospitals and Labs




MAGLUMI® X3

Features & Benefits

- Small but powerful, the throughput is up to 200 tests/hour, and throughput per unit area is 294 T/h/m².
- Compatible with all MAGLUMI® reagents, one of the broadest automated CLIA test menus in the world (166 parameters).
- The latest intelligent washing technology and bidirectional temperature control measurement guarantee accurate and reliable results. The comprehensive advanced design of MAGLUMI® X3 ensures excellent performance.
- Single reaction cup can avoid light pollution and increase cuvette utilization, its integrated packaging can avoid the stuck of the cuvette and destroying the cuvette wall.
- No-pause loading/unloading of reagents/samples without waiting or interrupting tests.

Functional Modules





All Balanced and Strong





Reagent and Sample Loading

- 72 sample positions with barcode reader
- Sample in each position can be defined as STAT
- 20 reagent positions with RFID reader, 24h refrigerated at 10°C (±0.3°C)
- No-pause loading/unloading of reagents/samples without waiting or interrupting tests
- Reagent and sample indicator light available, real-time status monitoring, no need to check on the computer monitor



Single Reaction Cup with Integrated Packaging

- Single reaction cup can avoid light pollution and increase the cuvette utilization
- Single reaction cup with integrated packaging can avoid the stuck of the cuvette and destroying the cuvette wall



Intelligent UI Interface

- New user-friendly UI interface (Support multi-language)
- Real-time display for monitoring each test procedure, reagent, and consumable quantities, it's convenient for the mangement
- Intelligent error recovery



Convenient Management of Analyzer Status

 Intuitive indicator light of reagent, sample, and consumables, no need to focus on the monitor, the status can be known by looking at the analyzer from a distance

Accurate Pipetting Technology

- Single needle coated with TEFLON to avoid carryover
- Independent washing unit for the needle, with different washing modes depending on different assay
- Pipetting detection technologies include liquid level, clot, and crash detection can ensure accurate and normal pipetting



Accurate Incubation Technology

- Single reaction cup with five-sided heating at 37.0°C (±0.3°C), ensuring accurate, quick and uniform incubation
- Incubating 80 cuvettes simultaneously
- Non-contact mixing with different modes before incubation, ensuring sufficient reaction




High-Efficiency Washing Technology

- Magnetic separation washing with 4-step independent units and non-contact vortex mixing can reduce non-specific binding
- Different mixing modes depending on different assays, ensuring a sufficient washing to improve sensitivity



Stable and Accurate Measuring Chamber

- ABEI-label exclusive proprietary technology
- PMT with heating and cooling function, ensuring the detection precision
- Independent measuring room, avoiding light cross-contamination



Outstanding Technologies



Exclusive Intelligent Washing Technology

Magnetic separation washing and non-contact vortex mixing






Non-contact vortex mixing units in the wash station, with different mixing modes depending on different assays, ensuring a sufficient washing to reduce the non-specific binding.


Intelligent washing machine




Exclusive intelligent washing technology can be understood as how the intelligent washing machine work. Different washing modes can be selected according to different clothes materials.

Stable and Accurate Measuring Chamber







Automatic temperature control in the measuring chamber can keep a stable background value, not affected by the ambient temperature, ensuring the accuracy of measurement results.



Top 10 reasons why choose MAGLUMI® X3

- One of the highest space efficiency CLIA analyzers combines one of the broadest automated CLIA test menus.
- The throughput is up to 200 tests/hour (Throughput per unit area is 294 T/h/m²). No worry about the lack of laboratory space.
- MAGLUMI® reagent kits (166 parameters) and consumables are compatible with MAGLUMI® X3.
- The comprehensive advanced design of MAGLUMI® X3 ensures excellent performance.
- Single reaction cup can avoid light pollution and increase cuvette utilization, its integrated packaging can avoid the stuck of the cuvette and destroying the cuvette wall.
- No-pause loading/unloading of reagents/samples without waiting or interrupting tests.
- Pipetting detection technologies include liquid level, clot, and crash detection can ensure accurate and normal pipetting.
- The latest intelligent washing technology and bidirectional temperature control measurement guarantee accurate and reliable results.
- Intelligent UI interface, user-friendly, easy, and convenient operation.
- Intuitive indicator light make users no need to check reagents and consumables frequently.

MAGLUMI® X3 Reagent and Consumables







Reaction Module



Starter 1+2



Light Check



– Wash Concentrate –



- System Tubing Cleaning Solution -

Broad CLIA Test Menu with 166 Assays



Thyroid

TSH (3rd Generation) T4

T3 FT4

FT3 Tg (Thyroglobulin) TGA (Anti-Tg) Anti-TPO TRAb TMA

Rev T3
*T-Uptake

Hepatic Fibrosis

HA
PIIIP N-P
C IV
Laminin
Cholyglycine

TORCH

Toxo IgG
Toxo IgM
Rubella IgG
Rubella IgM
CMV IgG
CMV IgM
HSV-1/2 IgG
HSV-2 IgG
*HSV-2 IgM
*HSV-1 IgG
*HSV-1 IgG
*HSV-1 IgG

EBV

EBV EA IgG EBV EA IgA EBV VCA IgG EBV VCA IgM EBV VCA IgA EBV NA IgG EBV NA IgA

Inflammation Monitoring

hs-CRP PCT (Procalcitonin) IL-6 (Interleukin 6) SAA(Serum Amyloid A)

Fertility

FSH LH HCG/β-HCG PRL (Prolactin) Estradiol Testosterone free Testosterone DHEA-S Progesterone free Estriol 17-OH Progesterone AMH SHBG Androstenedione *PIGF *sFlt-1

Hypertension

Direct Renin Aldosterone Angiotensin I Angiotensin II Cortisol ACTH

Autoimmune

Anti-CCP Anti-dsDNA IgG ANA Screen **ENA Screen** Anti-Sm IaG Anti-Rib-P IgG Anti-Scl-70 IgG Anti-Centromeres IgG Anti-Jo-1 IgG Anti-M2-3E IgG Anti-Histones IgG Anti-nRNP/Sm IgG Anti-SS-B IgG Anti-SS-A IgG TGA(Anti-Tg) Anti-TPO **TRAb** TMA ICA IAA(Anti Insulin) GAD 65 Anti-IA2 *Anti-Cardiolipin IgG *Anti-Cardiolipin IgM

Immunoglobulin

IgM IgA IgE IgG

*Anti-MPO

Tumor Markers

AFP

CEA Total PSA f-PSA CA 125 CA 15-3 CA 19-9 PAP CA 50 CYFRA 21-1 CA 242 CA 72-4 NSE S-100 SCCA TPA-snibe **ProGRP** HE4 HER-2 PIVKA-II

Infectious Disease

HBsAg Anti-HBs **HBeAg** Anti-HBe Anti-HBc Anti-HCV Syphilis Anti-HAV HAV IgM HIV Ab/Ag Combi Chagas HTLV I+II H.pylori IgG H.pylori IgA H.pylori IgM 2019-nCoV IgG 2019-nCoV IaM SARS-CoV-2 S-RBD IgG SARS-CoV-2 Neutralizing Antibody SARS-CoV-2 Ag *Anti-HBc IgM

Drug Monitoring

Digoxin CSA (Cyclosporine A) FK 506 (Tacrolimus)

Kidney Function

β₂-MG Albumin *NGAL

Cardiac

CK-MB Troponin I Myoglobin hs-cTnI H-FABP NT-proBNP BNP D-Dimer Lp-PLA2 *MPO

Metabolism

Pepsinogen I Pepsinogen II Gastrin-17 GH (hGH) IGF-I IGFBP-3

Prenatal Screening

AFP(Prenatal Screening)
Free β-HCG
PAPP-A
free Estriol

Anemia

Vitamin B12 Ferritin Folate (FA) *RBC Folate *FPO

C-Peptide

Glyco Metabolism

Insulin GAD 65 Anti-IA2 ICA IAA (Anti Insulin) Proinsulin

Bone Metabolism

Calcitonin
Osteocalcin
25-OH Vitamin D
Intact PTH
*β-CrossLaps (β-CTx)
*total P1NP

* Available soon

MAGLUMI® X3

Fully Automatic Chemiluminescence Immunoassay (CLIA) System

Technical Specifications

Chemiluminescence Principle

- Direct chemiluminescence
- Magnetic microbeads separation, fast and efficient
- Key technology of ABEI-label and magnetic microbead separating enhances stability and sensitivity of MAGLUMI® reagents

Main Features

- Throughput: up to 200 tests/hour
- Space occupied<0.68m², throughput of per unit area>294T/h/m²
- 24 hours ready to use

Reagent Loading

- 20 reagent positions on board with no-pause loading/unloading
- 24-hour refrigerated, constant at 10°C (±0.3°C) in reagent area
- · RFID reading all information of reagent
- Indicator light available, real-time status monitoring

Sample Loading

- Up to 72 sample positions
- No-pause loading/unloading, STAT available
- Indicator light available, real-time status monitoring
- Barcode recognition: Barcode reader available
- Types of sample tube: micro sample cup, blood collection tube, plastic test tube
- Sample editing mode:Lis connection, barcode recognition, manual editing

Reagent Features

- MAGLUMI® reagent kits are compatible with all MAGLUMI® analyzers
- Calibrators and internal controls included
- 10 points master curve, 2-points re-calibration
- RFID tag storing all informations of reagents
- Storage temperature 2°C-8°C

Reaction System

- Single reaction cup with five-sided heating at 37.0°C (±0.3°C) incubation, ensuring accurate, quick and uniform incubation
- Incubating 80 cuvettes simultaneously
- Non-contact mixing with different modes before incubation, ensuring sufficient reaction

 $\label{eq:maglumi} \mathbf{MAGLUMI^{\$}}, \mathbf{Snibelis^{TM}} \ \text{and} \ \mathbf{SnibeLinker^{TM}} \ \text{are trademarks of Snibe}.$ All other product names and trademarks are the property of their respective owners.

Sampling System

- Titanium needle for crushproof
- Teflon coated for carry-over prevention
- Liquid level detection, clot detection, prevent bubble interference, collision detection, wash automatically

Consumable System

- Single cup design, no-pause loading/unloading
- Starters: 230mL of Starter 1 and Starter 2
- System liquid: Up to load 10L at one time
- Indicator light monitors real-time consumable quantity

Washer

- 4-step washing technology
- Intelligent magnetic field design ensure excellent separation
- Three independent non-contact magnetic microbeads vortex cleaning technology
- Resuspension washing mode efficiently reduces non-specific binding and improves sensitivity
- Independent washing unit for waste liquid needle can completely wash without residue

Measuring Chamber

- High-sensitivity and low-noise photomultiplier tube (PMT)
- Heating and cooling bidirectional temperature control function included
- Original light pollution isolation technology prevents light pollution and ensures accurate and reliable results

Software System

- Friendly software UI design, humanized instrument operation procedure
- Real-time status display for each test procedure
- Self-recovery function make sure the analyzer running smoothly without unnecessary manual intervention
- \bullet Bidirectional $\;Snibelis^{\scriptscriptstyle TM}\;$ connection by TCP/IP and COM Port
- Supporting SnibeLinkerTM remote diagnosis system

Dimensions & Weight

- Dimensions: 78 (h) \times 90 (w) \times 75 (d) cm
- Weight: 127 kg









Shenzhen New Industries Biomedical Engineering Co., Ltd.