

CRITICOOL[®] Safe, Precise, and Effective^{3,2}

PATIENT TEMPERATURE MANAGEMENT



CritiCool CritiCool

Q

— NEXT GenWorks — Defining tomorrow, *today*.

TARGETED **TEMPERATURE** MANAGEMENT







CRITICOOL®

CritiCool[®] is a thermal regulating system indicated for maintaining pre-set body temperature as determined by a physician. The system is composed of two elements, the CritiCool device and the CureWrap[™].

The CritiCool device functions as a control unit for use with the CureWrap. The control unit reads the patient's core temperature through specific probes every 66 milliseconds and adjusts temperature as necessary in order to maintain the physician determined pre-set temperature settings.

The cooling/heating pump brings the water to the temperature set point and circulates it in a closed-loop system through the uniquely designed CureWrap.

FAST. SIMPLE. PRECISE.

- Accurate temperature control as set by the physician
- Full-color, touch-screen display
- Simple, easy to use interface with intuitive controls
- Real-time monitoring and graphical display of temperature measurements
- Quick and easy setup
- Track data using the CliniLogger[™]

CUREWRAP[™]


The CureWrap is a single use, one-piece flexible garment that acts in conjunction with the CritiCool control unit. The CureWrap's flexible configurations facilitate the wrapping of individual parts of the body to maximize surface coverage.


The CureWrap is designed to conveniently attach to the control unit using Quick Coupling Connectors (QCC) that allow for inflow and outflow. Integrated channels within the suit-like garment circulate water in a closed-loop system over a large surface area.

The CureWrap is secured via pressure sensitive adhesive strips which adhere to the garment, and not the patient.

DESIGNED WITH YOUR PATIENT IN MIND

- One-piece body shaped garment with various sizes and configurations
- Up to 85% body surface area coverage
- Soft mesh that is easy to wrap and secure
- Clutter-free with only one pair of connecting tubes







Soft-mesh garment with

built-in pressure relief channels



WRAP TYPE

Infant


Pediatric (Small) Pediatric (Medium) Pediatric

Pediatric (X-Large)

(Large)

Adult

TARGETED TEMPERATURE MANAGEMENT HAS NEVER BEEN EASIER



The CureWrap is available in a range of sizes and designs to accommodate your patient's anatomy:

РАСК ТҮРЕ	PART NUMBER	QUANTITY PER BOX	PATIENT HEIGHT/ WEIGHT	WRAP LENGTH/ WIDTH (M)
Single Size	508-03518	8	2.5-4 kg (5.5- 8.8 lb)	0.659 / 0.448
Single Size	508-03521	8	4-7 kg (8.8-15.4 lb)	0.698 / 0.602
	PED-SM008	8 (multi-size)	(x4) 2.5-4 kg (5.5- 8.8 lb) (x4) 4-7 kg (8.8-15.4 lb)	0.659 / 0.448 0.698 / 0.602
Assorted	PED-MD008	8 (multi-size)	(x4) 7–11 kg (15.4–24.3 lb) (x4) 79–91 cm (31–36")	0.981 / 0.628 1.118 / 0.740
Sizes	PED-LA008	8 (multi-size)	(x4) 91-104 cm (36-41") (x4) 104-122 cm (41-48")	1.225 / 0.841 1.390 / 1.054
	PED-XL008	8 (multi-size)	(x4) 122-135 cm (48-53") (x4) Over 135 cm (Over 53")	1.582 / 1.193 2.030 / 1.354
Single Size	508-03500	8	Over 135 cm (Over 53") Fits most adults	2.030 / 1.354

Stafford et al.: Therapeutic Hypothermia During Neonatal Transport: Active Cooling Helps Reach the Target. THERAPEUTIC HYPOTHERMIA AND TEMPERATURE MANAGEMENT 2016. ver et al.: Intravascular versus surface cooling for targeted temperature management after out-of-hospital cardiac arrest – an analysis of the TTM trial data. Critical Care 2016 20:381

ORDERING INFORMATION

CONTROL UNIT	
CritiCool® (115V)	Part #100-00003
CritiCool® (230V)	Part #100–00004
CritiCool® (100V)	Part #100-00007
CRITICOOL ACCESSORY KIT (Select o	ne only for each CritiCool
INFANT KIT FOR DISPOSABLE	Part #200-00330
TEMPERATURE SENSORS	
The Infant Kit contains:	
Core Adapter Cable	Part #014-00028
 Surface Adapter Cable 	Part #014-00129
 Connecting Water Tubes (2 by 2 way) 	Part #200-00109
ADULT KIT FOR DISPOSABLE	Part #200-00310
TEMPERATURE SENSORS	
The Adult Kit contains:	
Core Adapter Cable	Part #014-00028
 Surface Adapter Cable 	Part #014-00129

• Connecting Water Tubes (2 by 3 way)

Note: Disposable Temperature Probes must be ordered separately

REUSABLE KITS: (NOT AVAILABLE IN THE UNITED STATES AND SELECT MARKETS)

INFANT REUSABLE KIT

- The Infant Kit contains:
- Reusable Core Temperature Probe, Infant Part #014-00005
- Reusable Surface Temperature Probe
- Filter unit
- Connecting Water Tubes (2 by 2 way)

ADULT REUSABLE KIT

- The Adult Kit contains:
- Reusable Core Temperature Probe, Adult Part #014-00020
- Reusable Surface Temperature Probe
- Filter unit
- Connecting Water Tubes (2 by 3 way)

ACCESSORIES

Disposable Core Temperature Probe Disposable Surface Temperature Probe CliniLogger™

Note: Connecting Water Tubes (2 by 3 way) (Part #200-00147) is required for Pediatric X-Large CureWrap (Part #PED-XL008) and for Adult CureWrap (Part #508-03500).

HARDWARE SPECIFICATIONS

HEAT EXCHANGERS

Peltier Technology (TECs)

WATER TANK

Sterile or 0.22 micron filtered tap water Tank Capacity: 6 liters (1.6 gal.)

WATER TEMPERATURE

Water Temperature Accuracy: ± 0.3°C (0.54°F) Water Temperature (outflow): 13-40.8°C (55.4-105.4°F)

PUMP

unit)

Water Circulating Pump Pump Rate: 1.2 L/min Protected by 263 micron filter

PATIENT TEMPERATURE CHANNELS

2 Channels: Core, Surface Sensor Temperature Range: 15°C to 44°C (59°F to 111.2°F) Sensor Temperature Accuracy: ± 0.3°C (0.45°F)

TEMPERATURE AND PRESSURE SENSORS

System Sensors:

- 3 Internal Temperature Sensors: Water in, Water out, Thermostat
- 2 Pressure sensors
- Safety Measures:
- Over pressure protection and alarm
- High water temperature protection and alarm

PHYSICAL DIMENSIONS

Mobile unit with 4 wheels and 2 brakes 260 mm W x 625 mm D x 940 mm H / (10.23"W x 24.6"D x 37"H)

NET WEIGHT

34 kg / 75 lb

INPUT POWER

115VAC 5.7 A 50/60 HZ 230VAC 2.9 A 50/60 HZ 100VAC 6.6 A 50/60 HZ

LCD DISPLAY

Size: 5.7" color display Resolution: 320×240

CUREWRAP™

Single Patient Use

For more information about the products we offer, contact us at sales@belmontmedtech.com.









MANAGEMENT



Contact Information: Koteswaran S Genworks Health Pvt Ltd Unit 522-524, 5th Floor, Gamma Block Sigma Tech Park, Ramagondanahalli Varthur Road, Bangalore - 560066 Ph no. +91 9787143194 koteswaran.s@genworkshealth.com

Belmont Instrument LLC Certified to ISO 13485:2016 File Number: A18173

Belmont Medical Technologies[™] is a trademark of Belmont Instrument, LLC 2019 © Belmont Instrument, LLC | 701–00273 Rev. A





Part #014-00322 Part #014-00321 Part #017-00250

Part #200-00300

Part #200-00147

Part #200-00320

Part #014-00021

Part #200-R0130

Part #200-00109

Part #200-00147

Part #014-00021 Part #200-R0130