

## ONE SYSTEM. YOUR COMPLETE SOLUTION.

Together, the EleVision™ HD 1 system components enable — and empower surgeons to optimize minimally invasive surgeries (MIS)

#### 50 L insufflator ———

- Simple-to-use pressure and flow settings<sup>2</sup>
- Easy-to-read display of actual pressure and flow values<sup>2</sup>
- Individual programming options<sup>2</sup>
- Offers gas heating to target a 37 C body temperature<sup>3</sup>

#### CMOS technology camera -

- Captures native, HD images
- Maintains sharpness
- Provides better resolution<sup>4</sup>
- Enables color enhancement that can improve identification and characterization of critical structures<sup>1</sup>

#### Laparoscopic suction irrigation pump

- Radio Frequency Identification (RFID) control
- 21/m suction irrigation



Monitor
Displays 1080p
HD video



#### Laparoscope

Features a laserwelded design that lasts for more than 2,000 autoclave cycles<sup>5,†</sup>

### **LED** light source

- Features 20,000hour lifespan<sup>6</sup>
- Delivers 40 times more life than a Xenon<sup>™</sup> light bulb<sup>6</sup>

**Medtronic** 

Further, Together





COMPANY	MODEL	IMAGE ENHANCEMENT MODE(S)	CMOS/ CCD	LIGHT SOURCE
Medtronic	EleVision™ HD 1 System	Smoke reduction Color enhancement	CMOS	LED
Karl Storz™*	Image 1 S Connect II (CCU) Image 1 H3-Z FI (camera head)	Chroma Clara	CCD	LED
Olympus™*	CV-190 (CCU) CH-S190-XZ-E/Q (camera head)	NBI	CCD	Xenon <sup>™*</sup>
Stryker™*	1488	Clarity (not included) InfraVision (not included)	CMOS	LED (L9000)
	1288	Clarity (not included) InfraVision (not included)	CCD	Xenon <sup>™*</sup> (X8000)
	1188	None	CCD	Xenon™* (X8000)
Richard Wolf™*	3-Chip HD camera head Endocam Logic	Color contrast I-III	CCD	LED
	1-Chip HD camera head Endocam Logic	Color contrast I-III	CMOS	LED
ConMed™*	True HD 3CMOS CCU True HD 3CMOS camera head	None	CMOS	LED

# INSUFFLATOR COMPETITIVE COMPARISON

COMPANY	MAX GAS FLOW L/MIN	GAS WARMING	FEATURES
Medtronic 50 L Insufflator	50	YES	Pediatric, vessel harvesting, and bariatric modes for precise insufflation
Storz™* Endoflator 50	50	YES	Sensitive mode for sensitive applications (e.g., in pediatric surgery)
Olympus™* UHI-4	45	NO	Automatic smoke evacuation
Richard Wolf™* Highflow 45	45	YES	Four different hardware versions
Stryker™* Pneumosure	45	YES	Real-time pressure sensing
Aesculap™* Flow 40	45	YES	Gas warming

†Based on competitive website Aug 13, 2019 and competitive IFU.



#### Genworks Health Pvt Ltd

Unit 522-524, 5th Floor, Gamma Block, Sigma Tech Park, Ramagondanahalli, Varthur Road, Bangalore – 560066

Contact Info:-

Pooja Uday Parulekar

©+91 96422 14646

☑pooja.parulekar@genworkshealth.com







