



SOLVING CERVICAL ISSUES HAS NEVER BEEN EASIER

High-quality screening & High-precision treatment Screen & Treat in a single sitting!

Elimination of cervical cancer by 2030



In 2019, WHO issued new guidelines for the use of thermal ablation for cervical pre-cancerous lesions, recommending the thermal ablation as an alternative method of treatment.

















Treat up to CIN 1 & 2 lesions. Cervicitis. Ectropion/Ectopy or any other cervical abnormality.

ThermoGlide

An immediate point-of-care treatment option to prevent cervical issues

A complete thermo coagulation solution for the treatment of precancerous cervical lesions

The ThermoGlide uses heat to decimate tissue. The probe is heated quickly within **8 seconds and remains at 100°C** allowing for rapid cervical ablation. The superficial epithelium withers and dissipates after the treatment, while the stroma and glandular crypts are destroyed.

Thermo coagulation offers shorter treatment time as compared to cryotherapy.



- 4 probes
 Versatile treatment options
- ② 2 battery packs
 Wireless battery pack
- Charger

 Rechargeable battery, with no need for electricity
- Minimal pain and discomfort
 Compared to the conventional treatment
 methods
- Timer function & constant temperature setting reducing possibility of human error
- Quick learning curve, simple handling and assembly



In 2019, WHO issued new guidelines for the use of thermal ablation for cervical pre-cancerous lesions, recommending the thermal ablation as an alternative method of treatment.

Pre-cancerous lesions identified by EVA PRO can now be treated on site